

Gradient-Index Research Group Metrology Instruments

Instrument	Wavelength (μm)	Measurement		Material Type	Comments	Sample Geometry	Sample Prep Requirements
		Type	Index Range				
Visible Mach-Zehnder Interferometer	0.35 – 1.1	Relative Index	Any	GRIN	Gives two-dimensional relative index. Accuracy limited by sample preparation and mechanical properties.	Thin (~1mm), plane, parallel slice. Index is constant through sample thickness (exposed gradient).	Figure error < 1λ . Wedge error < $5\mu\text{m}$ over region of interest.
Schmidt Immersion	0.35 – 1.1	Absolute Index	< 1.7	GRIN	Same as Mach-Zehnder but uses index fluid that matches the sample at some point. Limited by available fluids.	Thin (~1mm), plane, parallel slice. Index is constant through sample thickness (exposed gradient).	Figure error < 1λ . Wedge error < $5\mu\text{m}$ over region of interest.
MWIR Mach-Zehnder Interferometer	1.1 – 4.6	Relative Index	Any	GRIN	Gives two-dimensional relative index. Accuracy limited by sample preparation and mechanical properties. Cannot use index fluid.	Thin (~2mm), plane, parallel slice. Index is constant through sample thickness (exposed gradient).	Figure error < 1λ . Wedge error < $5\mu\text{m}$ over region of interest.
LWIR Mach-Zehnder Interferometer	9.3 – 11.3	Relative Index	Any	GRIN	Gives two-dimensional relative index. Accuracy limited by sample preparation and mechanical properties. Cannot use index fluid.	Thin (~3mm), plane, parallel slice. Index is constant through sample thickness (exposed gradient).	Figure error < 1λ . Wedge error < $5\mu\text{m}$ over region of interest.
Abbe and Pulfrich Refractometers	0.4 – 1.1	Absolute Index (Surface)	< 1.7	Homogeneous	Requires index fluid and well-characterized reference prism. Sample must have lower index than prism and fluid.	One polished surface (recommended > 1 cm^2) and one polish/fine ground surface at 90° ; minimum to no bevel.	Figure error < 3λ
Metricon Refractometer	0.4 – 1.7	Absolute Index (Surface)	< 2.0	Homogeneous	Critical angle measurement with reference prism.	One polished surface (> 1 cm^2) required	Figure error < 3λ
Visible Sagnac Interferometer	0.35 – 1.1	Absolute Index (Bulk)	Any	Homogeneous	Accuracy limited by sample thickness measurement and sample preparation.	Plane parallel window. Diameter > 15mm Thickness > 5mm	Wedge < 5 arcmin Thicker samples reduce measurement uncertainty
MWIR Sagnac Interferometer	1.1 – 4.6	Absolute Index (Bulk)	Any	Homogeneous	Accuracy limited by sample thickness measurement and sample preparation.	Plane parallel window. Diameter > 20mm Thickness > 5mm	Wedge < 5 arcmin Thicker samples reduce measurement uncertainty

*Instrument under construction or upgrade, expected completion by Fall of 2016

6/23/2016

Gradient-Index Research Group Metrology Instruments

LWIR Sagnac Interferometer	9.3 – 11.3	Absolute Index (Bulk)	Any	Homogeneous	Accuracy limited by sample thickness measurement and sample preparation.	Plane parallel window. Diameter > 25mm Thickness > 5mm	Wedge < 5 arcmin Thicker samples reduce measurement uncertainty
*Low Coherence Distance Measuring Interferometer	VIS broadband	Absolute Thickness	N/A	GRIN or Homogeneous	Measures thickness map of plane parallel samples. Expected accuracy around $\pm 100\text{nm}$.	Nominally plane parallel sample with polished surfaces	Wedge < 5 arcmin
Beam Deflection	0.45-0.66	Relative Index	Any	GRIN	Calculate relative GRIN profile by measuring laser deflection through GRIN material.	One-dimensional gradient. Plane parallel window or cylinder geometry	Wedge < 2°
*Thermal Interferometer	0.46-0.633 *3.39	dn/dT CTE	Any (-40 to 50°C) (*-40 to 80°C)	GRIN or Homogeneous	Gives two-dimensional measurement of index and thickness change as a function of temperature.	Thin (~1mm), plane, parallel slice. Index is constant through sample thickness (exposed gradient).	Figure error < 1λ Wedge error < 3μm over region of interest. Reflective coating on one half of one surface
*Mini Mach-Zehnder Interferometer	0.6328 *0.5435	Relative Index	Any	GRIN	Microscope configuration for small samples or large fringe densities.	Thin (~1mm), plane, parallel slice. Index is constant through sample thickness (exposed gradient).	Figure error < 1λ. Wedge error < 3μm over region of interest.

Instrument	Best Case Error Bar	Advantages	Disadvantages
Mach-Zehnder Interferometer	2×10^{-5}	<ul style="list-style-type: none"> - Very versatile and robust measurement method - Schmidt immersion method provides absolute index 	<ul style="list-style-type: none"> - Sample prep problems can limit measurement accuracy - Sample prep is typically destructive
Abbe and Pulfrich Refractometers	1×10^{-5}	<ul style="list-style-type: none"> - Robust measurement method 	<ul style="list-style-type: none"> - Limited by reference prism and available index fluids - Surface measurement only
Sagnac Interferometer	1×10^{-4}	<ul style="list-style-type: none"> - Average index measurement - Used for any index or wavelength range - Uses plane parallel sample rather than prism 	<ul style="list-style-type: none"> - Limited by sample thickness measurement - Sensitive to scatter and surface shape
Beam Deflection	1×10^{-3}	<ul style="list-style-type: none"> - Nondestructive test for certain geometries - Can be implemented in-process 	<ul style="list-style-type: none"> - Limited accuracy - Requires constant index in propagation direction
Thermal Interferometer	CTE and dn/dT to 2 or 3 digits	<ul style="list-style-type: none"> - Gives CTE and dn/dT as a function of (x,y) - Useful for GRIN and homogeneous materials 	<ul style="list-style-type: none"> - Double pass (reduces effective fringe resolution) - Requires environmentally controlled chamber

*Instrument under construction or upgrade, expected completion by Fall of 2016

6/23/2016